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


Human CD3E&CD3D Heterodimer Protein (Cat. No. **CDD-H52Wa**) on SDS-PAGE under reducing (R) and non-reducing (NR) condition and the purity of the protein is greater than 95%. The purity of the protein is more than 85% and around 80-90 kDa verified by SEC-MALS.

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Biotinylated Human CD3E&CD3D Heterodimer Protein, His,Avitag&Tag Free Batch consistency


Immobilized Biotinylated Human CD3E&CD3D Heterodimer Protein, His,Avitag&Tag Free (Cat. No. **CDD-H82W6**) at 1 µg/mL (100 µL/well) on Streptavidin (Cat. No. **STN-N5116**) precoated (0.5 µg/well) plate, can bind Monoclonal Anti-Human CD3 Antibody, Mouse IgG2a (Clone: OKT3) (Cat. No. **CDE-M120a**) with a linear range of 0.2-6 ng/mL (QC tested), and batch differences EC50<0.0001 µg/mL.

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
Biotinylated Human CD3E&CD3D Heterodimer Protein, His,Avitag&Tag Free ELISA  
0.1 µg of Biotinylated Human CD3E&CD3D Heterodimer Protein, His,Avitag&Tag Free per well

Immobilized Biotinylated Human CD3E&CD3D Heterodimer Protein, His,Avitag&Tag Free (Cat. No. **CDD-H82W6**) at 1 µg/mL (100 µL/well) on streptavidin precoated (0.5 µg/well) plate, can bind Monoclonal Anti-Human CD3 Antibody, Mouse IgG2a (Clone: OKT3) (Cat. No. **CDE-M120a**) with a linear range of 0.2-6 ng/mL.


Human CD3E&CD3D Heterodimer Protein, His Tag&Tag Free ELISA  
0.1 µg of Human CD3E&CD3D Heterodimer Protein, His Tag&Tag Free per well


Immobilized Human CD3E&CD3D Heterodimer Protein, His Tag&Tag Free (Cat. No. **CDD-H52W1**) at 1 µg/mL (100 µL/well) can bind UCHT1 with a linear range of 0.2-6 ng/mL.

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Human CD3E&CD3D Heterodimer Protein, His Tag&Tag Free ELISA  
0.1 µg of Human CD3E&CD3D Heterodimer Protein, His Tag&Tag Free per well


Immobilized Human CD3E&CD3D Heterodimer Protein, His Tag&Tag Free (Cat. No. **CDD-H52W1**) at 1 µg/mL (100 µL/well) can bind Monoclonal Anti-Human CD3 Antibody, Mouse IgG2a (Clone: OKT3) (Cat. No. **CDE-M120a**) with a linear range of 0.2-6 ng/mL.

Biotinylated Human CD3E&CD3D Heterodimer Protein, His,Avitag&Tag Free (Cat. No. **CDD-H82W6**) captured on Biotin CAP - Series S sensor Chip can bind Anti-Human CD3 Antibody, Mouse IgG2a (Clone: OKT3) with an affinity constant of 22.5 nM as determined in a SPR assay (Biacore T200).

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Immobilized Human BCMA, His Tag (Cat. No. [BCA-H522y](#)) at 2 µg/mL, add increasing concentrations of Bispecific T cell Engager (CD3 X BCMA) in 10% human serum and then add Biotinylated Human CD3E&CD3D Heterodimer Protein, His,Avitag&Tag Free (Cat. No. [CDD-H82W6](#)) at 0.2 µg/mL. Detection was performed using HRP-conjugated streptavidin with sensitivity of 15 ng/mL.

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分子	品番	種別	製品説明	
CD3E & CD3D	CDD-H5255	Human	Human CD3 epsilon&CD3 delta Heterodimer Protein, Fc Tag &Fc Tag (MALS verified)	New
	CDD-H52W1	Human	Human CD3 epsilon & CD3 delta Heterodimer Protein, His Tag &Tag Free (MALS verified)	Hot
	CDD-H82W6	Human	Biotinylated Human CD3 epsilon&CD3 delta Heterodimer Protein, His,Avitag™&Tag Free (MALS verified)	Hot
	CDD-H52Wa	Human	Human CD3 epsilon&CD3 delta Heterodimer Protein, Fc,His Tag&Fc,Flag Tag (MALS verified)	Hot
	CDD-H5258	Human	Human CD3E&CD3D Heterodimer Protein, Llama Fc&Llama Fc, low endotoxin	
	CDD-H82F6	Human	Biotinylated Human CD3 epsilon&CD3 delta Heterodimer Protein, Fc,His,Avitag™&Fc,Flag,Avitag™ (MALS verified)	
	CDD-H82W0	Human	Biotinylated Human CD3E&CD3D Heterodimer Protein, Fc,His Tag&Fc,Flag Tag, ultra sensitivity (primary amine labeling) (MALS verified)	
	CDD-M5219	Mouse	Mouse CD3 epsilon&CD3 delta Heterodimer Protein, His Tag& Flag Tag (MALS verified)	
	CDD-M52W2	Mouse	Mouse CD3 epsilon&CD3 delta Heterodimer Protein, Fc,His Tag&Fc,Flag Tag (MALS verified)	
	CDD-M82W5	Mouse	Biotinylated Mouse CD3 epsilon&CD3 delta Heterodimer Protein, His,Avitag™&Tag Free	
	CDD-R52D8	Rat	Rat CD3 epsilon&CD3 delta Heterodimer Protein, His Tag& Flag Tag (MALS verified)	
	CDD-R5257	Rat	Rat CD3 epsilon&CD3 delta Heterodimer Protein, Llama Fc Tag &Llama Fc Tag (MALS verified)	
	CDD-C5259	Cynomolgus	Cynomolgus CD3E&CD3D Heterodimer Protein, Llama Fc&Llama Fc, low endotoxin	
	CDD-C52W9	Cynomolgus	Cynomolgus CD3 epsilon&CD3 delta Heterodimer Protein, Fc, His Tag&Fc,Flag Tag (MALS verified)	
	CDD-C82W6	Cynomolgus	Biotinylated Cynomolgus CD3E&CD3D Heterodimer Protein, His,Avitag™&Tag Free (MALS verified)	
	CDD-C52W4	Cynomolgus	Cynomolgus CD3 epsilon&CD3 delta Heterodimer Protein, His Tag&Flag Tag (MALS verified)	
CDD-R52D7	Rabbit	Rabbit CD3 epsilon&CD3 delta Heterodimer Protein, His Tag &Flag Tag (MALS verified)		

分子	品番	種別	製品説明
CD3E & CD3D	CDD-R5259	Rabbit	Rabbit CD3E&CD3D Heterodimer Protein, Llama Fc&Llama Fc, low endotoxin (MALS verified)
CD3E & CD3G	CDG-H5253	Human	Human CD3 epsilon&CD3 gamma Heterodimer Protein, Fc Tag &Fc Tag (MALS verified) 
	CDG-H52W9	Human	Human CD3 epsilon&CD3 gamma Heterodimer Protein, Fc, His Tag&Fc,Flag Tag (MALS verified) 
	CDG-H52W5	Human	Human CD3 epsilon&CD3 gamma Heterodimer Protein, His Tag&Tag Free (MALS verified)
	CDG-M52W2	Mouse	Mouse CD3 epsilon&CD3 gamma Heterodimer Protein,Fc,His Tag&Fc,Flag Tag (MALS verified)
	CDG-M52D2	Mouse	Mouse CD3 epsilon&CD3 gamma Heterodimer Protein, His Tag &Flag Tag (MALS verified)
	CDG-M58D2	Mouse	Biotinylated Mouse CD3 epsilon&CD3 gamma Heterodimer Protein, His,Avitag™&Flag Tag (MALS verified)
	CDG-R52D2	Rat	Rat CD3 epsilon&CD3 gamma Heterodimer Protein, His Tag &Flag Tag
	CDG-R82D4	Rat	Biotinylated Rat CD3 epsilon&CD3 gamma Heterodimer Protein, His,Avitag™&Flag Tag (MALS verified)
	CDG-C5257	Cynomolgus	Cynomolgus CD3E&CD3G Heterodimer Protein, Fc Tag&Fc Tag (MALS verified)
	CDG-C52W6	Cynomolgus	Cynomolgus CD3E&CD3G Heterodimer Protein, His Tag&Tag Free (MALS verified)
	CDG-C8219	Cynomolgus	Biotinylated Cynomolgus CD3E&CD3G Heterodimer Protein, His,Avitag™&Flag Tag (MALS verified)
	CDG-RB52D5	Rabbit	Rabbit CD3 epsilon&CD3 gamma Heterodimer Protein, His Tag &Flag Tag (MALS verified)
	CDG-R82D3	Rabbit	Biotinylated Rabbit CD3 epsilon&CD3 gamma Heterodimer Protein, His,Avitag™&Flag Tag (MALS verified)
CD3 epsilon	CDE-H5223	Human	Human CD3 epsilon Protein, His Tag (MALS verified) 
	CDE-H5256	Human	Human CD3 epsilon Protein, Fc,His Tag
	CDE-H82E1	Human	Biotinylated Human CD3 epsilon Protein, His,Avitag™
	CDE-H8223	Human	Biotinylated Human CD3 epsilon Protein, His Tag, ultra sensitivity (primary amine labeling)
	CDE-HR2H4	Human	HRP-Human CD3 epsilon Protein, His Tag
	CDN-H5253	Human	Human CD3 epsilon Protein, Mouse IgG2a Fc Tag (MALS verified)
	CDE-M5256	Mouse	Mouse CD3 epsilon Protein, Fc Tag
	CDE-C5226	Cynomolgus	Cynomolgus CD3 epsilon Protein, His Tag (MALS verified)
	CDE-C5254	Cynomolgus	Cynomolgus CD3 epsilon Protein, Fc,His Tag
CD3 delta	CDD-H5224	Human	Human CD3 delta Protein, His Tag
	CDD-C52H8	Cynomolgus	Cynomolgus CD3 delta Protein, His Tag
CD3 gamma	CDA-H5254	Human	Human CD3 gamma / CD3G Protein, Fc Tag (MALS verified)

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

Bispecific T-cell Engager (CD3XBCMA) immobilized on CM5 Chip can bind Human CD3E&CD3D Heterodimer Protein, His Tag&Tag Free (Cat. No. CDD-H52W1) with an affinity constant of 31.8 nM as determined in a SPR assay (Biacore T200).

Loaded Human CD3XBCMA scFv on AHC Biosensor via DMF Filed Human BCMA, Fc Tag (Cat. No. BC7-H5254), can bind Cynomolgus CD3E&CD3D Heterodimer Protein, His Tag&Flag Tag (Cat. No. CDD-C52W4) with an affinity constant of 8 nM as determined in BLI assay (ForteBio Octet Red96e).

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