


SL-51. Na_v1.8 Channel Modulators

Voltage-gated sodium channels are involved in signal transduction in electrically excitable tissues such as those found in muscles, the heart, and nerves. Deregulated function of voltage-gated sodium channels is observed in various pathological and disease conditions including pain, CNS, and cardiac disorders. Specific channel isoforms, Na_v1.7, Na_v1.8 and Na_v1.9 are predominantly expressed in sensory neurons thus providing a strong rationale as anti-pain targets [1].

Based on the analysis of several known Na_v1.8-blockers we have developed a pharmacophore model which

has been used to carry out an exhaustive search on a 20K+ set of natural product-like compounds. The highest scoring compounds were screened *in vitro* in the human Na_v1.8/β3 sodium channel cell line. The Na_v1.8 blocking effects were evaluated using the IonWorks™ Barracuda system. Blockage of the sodium channel was performed at 1 μM of a tested compound using a stimulus voltage pattern measuring peak current amplitudes for two test pulses TP1 (tonic inhibition) and TP2 (inactivated state inhibition).



Signature Library 51

Formats	Supplementary Information
80 compounds per plate 0.1 mg; 1 mg; 2 mg dry film/powder 0.1 μmol; 1 μmol DMSO solutions	SL#51_SL#51_Nav18_ICh.sdf Supplementary Information Tonic block of Na _v 1.8, at 1.0 μM – TP1 Inactivated state inhibition of Na _v 1.8 at 1.0 μM – TP2

References:

1. Channels (Austin). 2015 Nov-Dec; 9(6): 360–366; doi: 10.1080/19336950.2015.1079674

Contact us:

USA: +1 336 721 1617
Japan: +81-80-3401-9097
Europe/Global:

mparisi@asinex.com
sota@asinex.com
lsadovenko@asinex.com