

# Apollo

## バルク在庫品リスト

| カタログNo.  | 品名                                                           | Cas No.    | Qty   | Qty2  | Qty3 | 備考                                                                                  | 保管条件                          |
|----------|--------------------------------------------------------------|------------|-------|-------|------|-------------------------------------------------------------------------------------|-------------------------------|
| BI141962 | n-Propyl Gallate (Ph. Eur.) pure, pharma grade               | 121-79-9   | 1Kg   |       |      | Harmful                                                                             |                               |
| BIA0972  | TBE buffer (10X)                                             |            | 1L    | 5L    | 10L  | Irritant/Teratogenic                                                                |                               |
| BIA1125  | D(+)-Sucrose (USP-NF, BP, Ph. Eur., JP) pure, pharma grade   | 57-50-1    | 1Kg   | 5Kg   |      | Irritant/Moisture Sensitive/Store under Argon                                       |                               |
| BIA1371  | Sodium Chloride (USP, BP, Ph. Eur., JP) pure, pharma grade   | 7647-14-5  | 1Kg   | 5Kg   |      | Irritant                                                                            |                               |
| BIA4692  | Light Liquid Paraffin                                        | 8042-47-5  | 1L    | 2.5L  |      | Harmful                                                                             |                               |
| BIA4831  | D-Mannitol low endotoxin (Ph. Eur., BP, USP, JP)             | 69-65-8    | 1Kg   |       |      | Irritant                                                                            |                               |
| BIB6269  | Semicarbazide hydrochloride                                  | 563-41-7   | 1Kg   |       |      | Toxic                                                                               |                               |
| BIB6303  | Iminodiacetic acid                                           | 142-73-4   | 1kg   | 2.5kg |      | Irritant                                                                            |                               |
| BIG0709  | Glycine >99%                                                 | 56-40-6    | 1Kg   | 5Kg   |      | Irritant                                                                            |                               |
| BIG1007  | Glycerol monostearate, purified                              | 31566-31-1 | 1Kg   | 5Kg   |      | Irritant                                                                            | Store at -20°C                |
| BIK0452  | Kojic acid 99%                                               | 501-30-4   | 1Kg   | 2.5Kg | 5Kg  | Harmful/Carcinogenic                                                                |                               |
| BIM3011  | TC 100 Insect Medium                                         |            | 1L    |       |      | Irritant                                                                            |                               |
| BIT0501  | Tween 80                                                     | 9005-65-6  | 1L    |       |      | Irritant                                                                            |                               |
| BIT0723  | DL-Tyrosine                                                  | 556-03-6   | 1Kg   |       |      | Irritant                                                                            |                               |
| BIT0803  | Tween 20                                                     | 9005-64-5  | 1L    |       |      | Irritant                                                                            |                               |
| BIT4110  | Thiourea                                                     | 62-56-6    | 1Kg   |       |      | Carcinogenic/Teratogenic/Light Sensitive/Store under Argon                          |                               |
| BIU4110  | Urea, Molecular biology grade 99.5%+                         | 57-13-6    | 1Kg   | 5Kg   |      | Irritant                                                                            |                               |
| DE50N    | Deuterium oxide >99.9 Atom % D                               | 7789-20-0  | 1Kg   |       |      | Irritant/Hygroscopic/Moisture Sensitive/Store under Argon                           | Packaged in screw top bottles |
| DE50P    | Deuterium oxide >99.9 Atom % D                               | 7789-20-0  | 1L    |       |      | Irritant/Hygroscopic/Moisture Sensitive/Store under Argon                           | Packaged in screw top bottles |
| IN1042   | Aluminium(III) potassium disulphate dodecahydrate, ACS grade | 7784-24-9  | 2.5Kg | 25Kg  |      | Irritant                                                                            |                               |
| IN1093   | Ammonium thiocyanate 98%+                                    | 1762-95-4  | 1kg   | 5kg   |      | Harmful/Hygroscopic/Air Sensitive/Store under Argon                                 |                               |
| IN2065   | Iron(III) chloride hexahydrate, ACS grade                    | 10025-77-1 | 1kg   | 5kg   |      | Corrosive/Harmful/Hygroscopic/Store under Argon                                     |                               |
| IN2086   | Iron(II) sulphate heptahydrate, ACS grade                    | 7782-63-0  | 1Kg   | 5Kg   |      | Harmful/Irritant/Air Sensitive/Hygroscopic/Store under Argon                        |                               |
| IN2440   | Magnesium acetate tetrahydrate 99%                           | 16674-78-5 | 1Kg   |       |      | Irritant/Hygroscopic/Moisture Sensitive/Store under Argon                           |                               |
| IN2713   | Nickel(II) sulphate hexahydrate 98%                          | 10101-97-0 | 1Kg   |       |      | Toxic/Carcinogenic                                                                  |                               |
| IN2794   | Oxalic acid dihydrate, ACS grade                             | 6153-56-6  | 1Kg   | 5Kg   |      | Corrosive/Harmful/Hygroscopic                                                       |                               |
| IN2901   | Potassium sorbate                                            | 24634-61-5 | 1Kg   |       |      | Irritant                                                                            |                               |
| IN2904   | Potassium dihydrogen phosphate 99.9%                         | 7778-77-0  | 1Kg   | 5Kg   |      | Irritant/Hygroscopic/Store under Argon                                              |                               |
| IN3685   | Titanium, sponge 0.8cm & down 99.5%                          | 7440-32-6  | 1Kg   | 5Kg   |      | Pyrophoric/Air Sensitive/Moisture Sensitive/Store under Argon                       |                               |
| IN3874   | Zinc, 99+%, mossy                                            | 7440-66-6  | 1Kg   | 5Kg   |      | Irritant/Moisture Sensitive/Store under Argon                                       |                               |
| OR0038   | 2,6-Dichloroaniline                                          | 608-31-1   | 1kg   |       |      | Toxic/Harmful/Irritant                                                              |                               |
| OR10004  | 1-Methyl-1H-imidazole 99%                                    | 616-47-7   | 1L    | 2.5L  |      | Toxic/Corrosive/Harmful/Hygroscopic/Store under Argon                               |                               |
| OR10026  | 3-Methylbutan-2-ol                                           | 598-75-4   | 1Kg   |       |      | Flammable/Harmful/Irritant                                                          |                               |
| OR10036  | 2-Methoxyphenol                                              | 90-05-1    | 1Kg   | 5Kg   |      | Harmful/Irritant/Light Sensitive/Air Sensitive/Moisture Sensitive/Store under Argon |                               |
| OR10815  | Phenyl dichlorophosphate 97%                                 | 770-12-7   | 1Kg   |       |      | Corrosive/Harmful/Moisture Sensitive/Store under Argon                              |                               |
| OR12500  | N-Methylbenzylamine 99%                                      | 103-67-3   | 1L    |       |      | Corrosive/Toxic/Harmful/Air Sensitive/Store under Argon                             |                               |
| OR13036  | Hydrobromic acid, 48% aqueous solution                       | 10035-10-6 | 1Kg   | 5Kg   |      | Corrosive/Irritant/Air Sensitive/Light Sensitive/Store under Argon                  |                               |
| OR14328  | Diphenylcarbonate                                            | 102-09-0   | 1Kg   | 10Kg  |      | Harmful                                                                             |                               |

【お問い合わせ】 ナミキ商事株式会社 医薬化学品部  
 Tel: 050-5527-9830 (部署直通)  
 Tel: 03-3354-4065 (代表) Fax: 03-3354-3527  
 mail: info@namiki-s.co.jp https://www.namiki-s.co.jp



# Apollo

## バルク在庫品リスト

| カタログNo.  | 品名                                       | Cas No.    | Qty   | Qty2  | Qty3 | 備考                                                                                                                            | 保管条件                        |
|----------|------------------------------------------|------------|-------|-------|------|-------------------------------------------------------------------------------------------------------------------------------|-----------------------------|
| OR16495  | 2-Ethylhexanoic acid                     | 149-57-5   | 2.5L  | 10L   |      | Harmful/Irritant                                                                                                              |                             |
| OR16612  | Triethyl orthoformate 98%                | 122-51-0   | 1L    | 2.5L  |      | Flammable/Irritant/Moisture Sensitive/Store under Argon/Keep Cold                                                             |                             |
| OR16937  | Allyl bromide 99%                        | 106-95-6   | 1kg   |       |      | Highly Flammable/Toxic/Corrosive/Carcinogenic/Mutagenic/Stench/Light Sensitive/Moisture Sensitive/Keep Cold/Store under Argon |                             |
| OR17127  | 1,4-Dibromobenzene                       | 106-37-6   | 1Kg   | 2.5Kg |      | Irritant/Light Sensitive/Store under Argon                                                                                    |                             |
| OR1880   | 2-Amino-1,3-thiazole 97%                 | 96-50-4    | 1kg   |       |      | Harmful/Irritant/Mutagenic/Keep Cold                                                                                          |                             |
| OR2121   | 3-Bromoaniline 98%                       | 591-19-5   | 1kg   |       |      | Toxic/Harmful/Irritant/Light Sensitive/Keep Cold                                                                              |                             |
| OR2477   | L-Threonine                              | 72-19-5    | 1kg   | 5kg   |      | Irritant                                                                                                                      |                             |
| OR25776  | Methyl(triphenyl)phosphonium bromide 99% | 1779-49-3  | 1Kg   |       |      | Toxic/Harmful/Hygroscopic/Store under Argon                                                                                   |                             |
| OR28814  | Ethyl formate                            | 109-94-4   | 5000g |       |      | Highly Flammable/Harmful/Irritant/Moisture Sensitive/Keep Cold/Store under Argon                                              |                             |
| OR40577  | 2-Aminobenzenesulphonic acid             | 88-21-1    | 1Kg   |       |      | Corrosive/Light Sensitive/Keep Cold                                                                                           |                             |
| OR4335   | 2-Bromobutane 98%                        | 78-76-2    | 1Kg   |       |      | Highly Flammable/Light Sensitive                                                                                              |                             |
| OR460052 | Tetrachloroethylene 99%                  | 127-18-4   | 1Kg   |       |      | Toxic/Carcinogenic                                                                                                            |                             |
| OR460066 | Citric acid monohydrate                  | 5949-29-1  | 1Kg   | 2.5Kg |      | Irritant/Hygroscopic/Store under Argon                                                                                        |                             |
| OR46222  | 6-Aminohexanoic acid                     | 60-32-2    | 1Kg   |       |      | Irritant                                                                                                                      |                             |
| OR52392  | Vanillin 97+%                            | 121-33-5   | 1kg   |       |      | Irritant/Light Sensitive/Air Sensitive/Moisture Sensitive/Store under Argon                                                   |                             |
| OR52524  | Benzyl cinnamate                         | 103-41-3   | 1kg   | 5kg   | 10kg | Irritant                                                                                                                      |                             |
| OR5344   | 2-Methoxybenzoic acid 98%                | 579-75-9   | 1Kg   | 5Kg   | 10Kg | Irritant                                                                                                                      |                             |
| OR5393   | Naphthalene-1,4-dione 97%                | 130-15-4   | 2.5kg | 5Kg   |      | Very Toxic/Toxic/Light Sensitive                                                                                              |                             |
| OR54367  | Methacrylic acid 99+%                    | 79-41-4    | 1Kg   | 5Kg   |      | Corrosive/Toxic/Harmful/Irritant/STENCH/Light Sensitive/Keep Cold/Store under Argon                                           | Stabilised with 250ppm MEHQ |
| OR54474  | Methyl acetoacetate                      | 105-45-3   | 1kg   | 5kg   |      | Irritant/Light Sensitive/Keep Cold                                                                                            |                             |
| OR5451   | Benzil                                   | 134-81-6   | 2.5Kg |       |      | Irritant                                                                                                                      |                             |
| OR54559  | Phthalimide potassium salt 99%           | 1074-82-4  | 1Kg   | 5Kg   |      | Irritant/Moisture Sensitive/Store under Argon                                                                                 |                             |
| OR5460   | L-(+)-Mandelic acid                      | 17199-29-0 | 1Kg   |       |      | Irritant/Corrosive/Moisture Sensitive/Light Sensitive/Store under Argon                                                       |                             |
| OR55016  | Lauric acid 98%                          | 143-07-7   | 1Kg   | 5Kg   | 25Kg | Irritant/Corrosive                                                                                                            |                             |
| OR55226  | Dimethyl adipate                         | 627-93-0   | 1Kg   |       |      | Irritant                                                                                                                      |                             |
| OR59808  | 2,2'-Oxydiethanol 99%                    | 111-46-6   | 2.5Kg | 10Kg  |      | Harmful/Hygroscopic/Store under Argon                                                                                         |                             |
| OR61047  | 4-Methylpentan-2-one 99%                 | 108-10-1   | 1Kg   |       |      | Highly Flammable/Harmful/Irritant                                                                                             |                             |
| OR61424  | 1,4-Dioxane 99%                          | 123-91-1   | 1L    | 2.5L  |      | Highly Flammable/Irritant/Carcinogenic/Hygroscopic/Air Sensitive/Store under Argon                                            |                             |
| OR62004  | Hydroxylamine hydrochloride 99%          | 5470-11-1  | 1Kg   |       |      | Corrosive/Toxic/Harmful/Irritant/Carcinogenic/Air Sensitive/Hygroscopic/Store under Argon                                     |                             |
| OR6566   | Pent-2-enoic acid                        | 626-98-2   | 1Kg   |       |      | Corrosive/Harmful                                                                                                             |                             |
| OR8791   | Chlorobenzene                            | 108-90-7   | 1L    |       |      | Flammable/Harmful                                                                                                             |                             |
| OR923775 | Calcium pantothenate                     | 137-08-6   | 1kg   |       |      | Irritant                                                                                                                      |                             |
| OR923999 | Benzyl alcohol                           | 100-51-6   | 1kg   | 10kg  |      | Irritant/Hygroscopic                                                                                                          |                             |
| OR9595   | 1-Bromooctane 98%                        | 111-83-1   | 1Kg   |       |      | Irritant                                                                                                                      |                             |
| PC1190   | Benzotrifluoride 99%                     | 98-08-8    | 1kg   |       |      | Highly Flammable                                                                                                              |                             |

【お問い合わせ】 ナミキ商事株式会社 医薬化学品部

Tel: 050-5527-9830 (部署直通)

Tel: 03-3354-4065 (代表) Fax: 03-3354-3527

mail: info@namiki-s.co.jp <https://www.namiki-s.co.jp>



# Apollo バルク在庫品リスト

| カタログNo. | 品名                                                    | Cas No.     | Qty     | Qty2 | Qty3 | 備考                                                                           | 保管条件           |
|---------|-------------------------------------------------------|-------------|---------|------|------|------------------------------------------------------------------------------|----------------|
| PC1210  | 3-Benzoyl-1,1,1-trifluoroacetone 98%                  | 326-06-7    | 1kg     |      |      | Irritant                                                                     |                |
| PC4382  | Fomblin® Y L-VAC 14/6                                 | 69991-67-9  | 1Kg     |      |      | Irritant                                                                     |                |
| PC4797  | Hexafluorosilicic acid, 23% aqueous solution          | 16961-83-4  | 1Kg     |      |      | Corrosive                                                                    |                |
| PC49140 | Methyl trifluoroacetate 97%                           | 431-47-0    | 1Kg     |      |      | Highly Flammable/Corrosive/Lachrymatory/Moisture Sensitive/Store under Argon |                |
| PC6294  | Potassium hexafluorosilicate 98%                      | 16871-90-2  | 1Kg     | 5Kg  |      | Toxic/Moisture Sensitive/Store under Argon                                   |                |
| PC6473  | Novac HFE-72DE Engineered Fluid (3M)                  | 839682-75-6 | 1Kg     | 5Kg  |      | Harmful/Irritant                                                             | Blend          |
| PMA267A | Anderson's basal salt mixture                         |             | 10 x 1L |      |      | Irritant                                                                     | Store at 2-8°C |
| PMA267B | Anderson's basal salt mixture                         |             | 10L     |      |      | Irritant                                                                     | Store at 2-8°C |
| PMA267C | Anderson's basal salt mixture                         |             | 50L     |      |      | Irritant/Moisture Sensitive/Store under Argon                                | Store at 2-8°C |
| PMB141A | BM-1 orchid medium                                    |             | 10 x 1L |      |      | Irritant/Moisture Sensitive/Store under Argon                                | Store at 2-8°C |
| PMB141C | BM-1 orchid medium                                    |             | 50L     |      |      | Irritant                                                                     | Store at 2-8°C |
| PMB141N | BM-1 orchid medium                                    |             | 10L     |      |      | Irritant/Keep Cold                                                           | Store at 2-8°C |
| PMB142A | BM-2 orchid medium                                    |             | 10 x 1L |      |      | Irritant                                                                     | Store at 2-8°C |
| PMB142C | BM-2 orchid medium                                    |             | 50L     |      |      | Irritant                                                                     | Store at 2-8°C |
| PMB142N | BM-2 orchid medium                                    |             | 10L     |      |      | Irritant                                                                     | Store at 2-8°C |
| PMC167A | Chu N6 basal salt medium with vitamins                |             | 10 x 1L |      |      | Oxidising Agent/Irritant/Hygroscopic/Store under Argon                       | Store at 2-8°C |
| PMC167B | Chu N6 basal salt medium with vitamins                |             | 10L     |      |      | Oxidising Agent/Irritant/Hygroscopic/Store under Argon                       | Store at 2-8°C |
| PMC167C | Chu N6 basal salt medium with vitamins                |             | 50L     |      |      | Oxidising Agent/Irritant/Hygroscopic/Store under Argon                       | Store at 2-8°C |
| PMC206A | Cape Sundew/Venus Flytrap Multiplication Basal Medium |             | 10 x 1L |      |      | Oxidising Agent/Irritant                                                     | Store at 2-8°C |
| PMC206B | Cape Sundew/Venus Flytrap Multiplication Basal Medium |             | 10L     |      |      | Oxidising Agent/Irritant                                                     | Store at 2-8°C |
| PMC206C | Cape Sundew/Venus Flytrap Multiplication Basal Medium |             | 50L     |      |      | Oxidising Agent/Irritant                                                     | Store at 2-8°C |
| PMC287  | Chee and Pool vitis Cd2 medium                        |             |         | 10L  |      |                                                                              |                |
| PMC416A | Chu N6 basal salt mixture                             |             | 10 x 1L |      |      | Irritant                                                                     | Store at 2-8°C |
| PMC416B | Chu N6 basal salt mixture                             |             | 10L     |      |      | Irritant                                                                     | Store at 2-8°C |
| PMC416C | Chu N6 basal salt mixture                             |             | 50L     |      |      | Irritant                                                                     | Store at 2-8°C |
| PMF522A | Fast Terrestrial Orchid Medium                        |             | 10 x 1L |      |      | Irritant/Moisture Sensitive/Hygroscopic/Store under Argon                    | Store at 2-8°C |
| PMF522B | Fast Terrestrial Orchid Medium                        |             | 10L     |      |      | Irritant/Moisture Sensitive                                                  | Store at 2-8°C |
| PMF522C | Fast Terrestrial Orchid Medium                        |             | 50L     |      |      | Irritant/Moisture Sensitive                                                  | Store at 2-8°C |
| PMG398A | Gamborg's B-5 basal salt medium                       |             | 10 x 1L |      |      | Irritant                                                                     | Store at 2-8°C |
| PMG398B | Gamborg's B-5 basal salt medium                       |             | 10L     |      |      | Irritant                                                                     | Store at 2-8°C |
| PMG398C | Gamborg's B-5 basal salt medium                       |             | 50L     |      |      | Oxidising Agent/Irritant/Keep Cold                                           | Store at 2-8°C |
| PMG768A | Gamborg's basal salt mixture                          |             | 10 x 1L |      |      | Irritant                                                                     | Store at 2-8°C |
| PMG768B | Gamborg's basal salt mixture                          |             | 10L     |      |      | Irritant                                                                     | Store at 2-8°C |
| PMG768C | Gamborg's basal salt mixture                          |             | 50L     |      |      | Oxidising Agent/Irritant/Hygroscopic/Store under Argon/Keep Cold             | Store at 2-8°C |
| PMH393  | Heller's basal salt mixture                           |             |         | 10L  |      |                                                                              |                |
| PMI365A | Ichihashi new phalaenopsis (NP) medium                |             | 10 x 1L |      |      | Irritant                                                                     | Store at 2-8°C |
| PMI365B | Ichihashi new phalaenopsis (NP) medium                |             | 10L     |      |      | Irritant                                                                     | Store at 2-8°C |
| PMI365C | Ichihashi new phalaenopsis (NP) medium                |             | 50L     |      |      | Irritant                                                                     | Store at 2-8°C |
| PMK400A | Knudson C orchid medium                               |             | 10 x 1L |      |      | Irritant                                                                     | Store at 2-6°C |
| PMK400B | Knudson C orchid medium                               |             | 10L     |      |      | Irritant                                                                     | Store at 2-6°C |
| PMK400C | Knudson C orchid medium                               |             | 50L     |      |      | Irritant                                                                     | Store at 2-6°C |
| PMK413  | Kao and Michayluk basal salt mixture                  |             |         | 10L  |      |                                                                              |                |
| PML154A | Lloyd and McCown's woody plant basal salt mixture     |             | 10 x 1L |      |      | Irritant/Moisture Sensitive                                                  | Store at 2-8°C |
| PML154B | Lloyd and McCown's woody plant basal salt mixture     |             | 1 x 10L |      |      | Irritant/Moisture Sensitive                                                  | Store at 2-8°C |
| PML154C | Lloyd and McCown's woody plant basal salt mixture     |             | 1 x 50L |      |      | Irritant/Moisture Sensitive                                                  | Store at 2-8°C |
| PML444A | Lloyd and McCown's woody plant micronutrient mixture  |             | 10 x 1L |      |      | Irritant                                                                     | Store at 2-8°C |

【お問い合わせ】 ナミキ商事株式会社 医薬化学品部

Tel: 050-5527-9830 (部署直通)

Tel: 03-3354-4065 (代表) Fax: 03-3354-3527

mail: info@namiki-s.co.jp https://www.namiki-s.co.jp



# Apollo

## バルク在庫品リスト

| カタログNo. | 品名                                                                                       | Cas No. | Qty     | Qty2 | Qty3 | 備考                                            | 保管条件           |
|---------|------------------------------------------------------------------------------------------|---------|---------|------|------|-----------------------------------------------|----------------|
| PML444B | Lloyd and McCown's woody plant micronutrient mixture                                     |         | 10L     |      |      | Irritant                                      | Store at 2-8°C |
| PML444C | Lloyd and McCown's woody plant micronutrient mixture                                     |         | 50L     |      |      | Irritant                                      | Store at 2-8°C |
| PML449A | Lloyd and McCown's woody plant medium with vitamins                                      |         | 10 x 1L |      |      | Irritant                                      | Store at 2-8°C |
| PML449B | Lloyd and McCown's woody plant medium with vitamins                                      |         | 10L     |      |      | Irritant                                      | Store at 2-8°C |
| PML449C | Lloyd and McCown's woody plant medium with vitamins                                      |         | 50L     |      |      | Irritant                                      | Store at 2-8°C |
| PMM153  | Murashige and Skoog modified basal salts (½ micros and ½ macros)                         |         |         | 10L  |      | Oxidising Agent/Keep Cold                     |                |
| PMM290  | Murashige and Skoog modified basal salt mixture                                          |         |         | 10L  |      |                                               |                |
| PMM401A | Murashige and Skoog modified basal medium (with 6-BAP)                                   |         | 10 x 1L |      |      | Irritant                                      | Store at 2-8°C |
| PMM401B | Murashige and Skoog modified basal medium (with 6-BAP)                                   |         | 10L     |      |      | Irritant                                      | Store at 2-8°C |
| PMM401C | Murashige and Skoog modified basal medium (with 6-BAP)                                   |         | 50L     |      |      | Irritant                                      | Store at 2-8°C |
| PMM407A | Murashige and Skoog modified basal salts (contains no N, P or K)                         |         | 10 x 1L |      |      | Irritant                                      | Store at 2-8°C |
| PMM407B | Murashige and Skoog modified basal salts (contains no N, P or K)                         |         | 10L     |      |      | Irritant                                      | Store at 2-8°C |
| PMM407C | Murashige and Skoog modified basal salts (contains no N, P or K)                         |         | 50L     |      |      | Irritant                                      | Store at 2-8°C |
| PMM462A | Musa (Banana) Multiplication Medium                                                      |         | 10 x 1L |      |      | Irritant                                      | Store at 2-8°C |
| PMM462B | Musa (Banana) Multiplication Medium                                                      |         | 10L     |      |      | Irritant                                      | Store at 2-8°C |
| PMM462C | Musa (Banana) Multiplication Medium                                                      |         | 50L     |      |      | Irritant                                      | Store at 2-8°C |
| PMM482A | Malmgren Modified Terrestrial Orchid Medium (without sucrose, agar and pineapple powder) |         | 10 x 1L |      |      | Irritant/Moisture Sensitive/Store under Argon | Store at 2-8°C |
| PMM482B | Malmgren Modified Terrestrial Orchid Medium (without sucrose, agar and pineapple powder) |         | 10L     |      |      | Irritant/Moisture Sensitive/Store under Argon | Store at 2-8°C |
| PMM482C | Malmgren Modified Terrestrial Orchid Medium (without sucrose, agar and pineapple powder) |         | 50L     |      |      | Irritant/Moisture Sensitive/Store under Argon | Store at 2-8°C |
| PMM502A | Murashige and Skoog modified basal salts (macro-nutrient salt base-1x)                   |         | 10 x 1L |      |      | Irritant                                      | Store at 2-8°C |
| PMM502B | Murashige and Skoog modified basal salts (macro-nutrient salt base-1x)                   |         | 10L     |      |      | Irritant                                      | Store at 2-8°C |
| PMM502C | Murashige and Skoog modified basal salts (macro-nutrient salt base-1x)                   |         | 50L     |      |      | Irritant                                      | Store at 2-8°C |
| PMM507A | Murashige Cattleya Orchid Multiplication Medium                                          |         | 10 x 1L |      |      | Irritant                                      | Store at 2-8°C |
| PMM507B | Murashige Cattleya Orchid Multiplication Medium                                          |         | 10L     |      |      | Irritant                                      | Store at 2-8°C |
| PMM507C | Murashige Cattleya Orchid Multiplication Medium                                          |         | 50L     |      |      | Irritant                                      | Store at 2-8°C |
| PMM508A | Murashige Fern multiplication medium                                                     |         | 10 x 1L |      |      | Irritant                                      | Store at 2-6°C |
| PMM508B | Murashige Fern multiplication medium                                                     |         | 10L     |      |      | Irritant                                      | Store at 2-8°C |
| PMM508C | Murashige Fern multiplication medium                                                     |         | 50L     |      |      | Irritant                                      | Store at 2-8°C |
| PMM517a | Murashige African violet/Gloxinia multiplication medium                                  |         | 10x1L   |      |      | Irritant                                      |                |
| PMM517b | Murashige African violet/Gloxinia multiplication medium                                  |         | 10L     |      |      | Irritant                                      |                |
| PMM517c | Murashige African violet/Gloxinia multiplication medium                                  |         | 50L     |      |      | Irritant                                      |                |
| PMM519A | Murashige and Skoog basal medium with vitamins                                           |         | 10x1L   |      |      | Oxidising Agent/Irritant                      | Store at 2-8°C |
| PMM519B | Murashige and Skoog basal medium with vitamins                                           |         | 10L     |      |      | Oxidising Agent/Irritant                      | Store at 2-8°C |
| PMM519C | Murashige and Skoog basal medium with vitamins                                           |         | 50L     |      |      | Oxidising Agent/Irritant                      | Store at 2-8°C |
| PMM522  | Murashige and Skoog modified basal salts (micro-nutrient salt base-1x)                   |         |         | 10L  |      |                                               |                |
| PMM524A | Murashige and Skoog basal mixture 10 x 1L                                                |         | 10 x 1L |      |      | Oxidising Agent/Irritant                      | Store at 2-8°C |
| PMM524B | Murashige and Skoog basal mixture 10L                                                    |         | 10L     |      |      | Oxidising Agent/Irritant                      | Store at 2-8°C |
| PMM524C | Murashige and Skoog basal mixture 50L                                                    |         | 50L     |      |      | Oxidising Agent/Irritant                      | Store at 2-8°C |
| PMM524D | Murashige and Skoog basal mixture 100L                                                   |         | 100L    |      |      | Oxidising Agent/Irritant                      | Store at 2-8°C |
| PMM527  | Murashige shoot tip rooting medium                                                       |         |         | 10L  |      |                                               |                |
| PMM531A | Murashige and Skoog modified basal salts (no nitrogen)                                   |         | 10 x 1L |      |      | Irritant                                      | Store at 2-8°C |
| PMM531B | Murashige and Skoog modified basal salts (no nitrogen)                                   |         | 10L     |      |      | Irritant                                      | Store at 2-8°C |

【お問い合わせ】 ナミキ商事株式会社 医薬化学品部

Tel: 050-5527-9830 (部署直通)

Tel: 03-3354-4065 (代表) Fax: 03-3354-3527

mail: info@namiki-s.co.jp <https://www.namiki-s.co.jp>



# Apollo

## バルク在庫品リスト

| カタログNo.  | 品名                                                                                  | Cas No. | Qty     | Qty2 | Qty3 | 備考                                             | 保管条件           |
|----------|-------------------------------------------------------------------------------------|---------|---------|------|------|------------------------------------------------|----------------|
| PMM531C  | Murashige and Skoog modified basal salts (no nitrogen)                              |         | 50L     |      |      | Irritant                                       | Store at 2-8°C |
| PMM535   | Murashige and Skoog modified basal medium                                           |         |         | 10L  |      |                                                |                |
| PMM541   | Murashige and Skoog modified basal salts (without KH <sub>2</sub> PO <sub>4</sub> ) |         |         | 10L  |      |                                                |                |
| PMM5501A | Murashige and Skoog medium with Vitamins and Sucrose                                |         | 10 x 1L |      |      | Oxidising Agent/Irritant                       | Store at 2-8°C |
| PMM5501B | Murashige and Skoog medium with Vitamins and Sucrose                                |         | 10L     |      |      | Oxidising Agent/Irritant                       | Store at 2-8°C |
| PMM5501C | Murashige and Skoog medium with Vitamins and Sucrose                                |         | 50L     |      |      | Oxidising Agent/Irritant                       | Store at 2-8°C |
| PMM5501D | Murashige and Skoog medium with Vitamins and Sucrose                                |         | 100L    |      |      | Oxidising Agent/Irritant                       | Store at 2-8°C |
| PMM551A  | Malmgren modified terrestrial orchid medium                                         |         | 10x1L   |      |      | Irritant/Moisture Sensitive/Store under Argon  | Store at 2-8°C |
| PMM551B  | Malmgren modified terrestrial orchid medium                                         |         | 10L     |      |      | Irritant/Moisture Sensitive/Store under Argon  | Store at 2-8°C |
| PMM551C  | Malmgren modified terrestrial orchid medium                                         |         | 50L     |      |      | Irritant/Moisture Sensitive/Store under Argon  | Store at 2-8°C |
| PMM555   | Murashige and Skoog modified basal medium (with Kinetin)                            |         |         | 10L  |      |                                                |                |
| PMM561   | Murashige and Skoog modified basal salts (½ Nitrogen)                               |         |         | 10L  |      |                                                |                |
| PMM571A  | Murashige and Skoog modified basal salts (without NH <sub>4</sub> NO <sub>3</sub> ) |         | 10 x 1L |      |      | Oxidising Agent/Irritant                       | Store at 2-8°C |
| PMM571B  | Murashige and Skoog modified basal salts (without NH <sub>4</sub> NO <sub>3</sub> ) |         | 10L     |      |      | Oxidising Agent/Irritant                       | Store at 2-8°C |
| PMM571C  | Murashige and Skoog modified basal salts (without NH <sub>4</sub> NO <sub>3</sub> ) |         | 50L     |      |      | Oxidising Agent/Irritant                       | Store at 2-8°C |
| PMM701   | Murashige and Skoog modified basal medium (with 2iP)                                |         |         | 10L  |      |                                                |                |
| PMM702   | Murashige and Skoog modified basal medium (with 2iP)                                |         |         | 10L  |      |                                                |                |
| PMN479A  | NLN medium                                                                          |         | 10 x 1L |      |      | Oxidising Agent/Irritant                       | Store at 2-8°C |
| PMN479B  | NLN medium                                                                          |         | 10L     |      |      | Oxidising Agent/Irritant                       | Store at 2-8°C |
| PMN479C  | NLN medium                                                                          |         | 50L     |      |      | Oxidising Agent/Irritant                       | Store at 2-8°C |
| PMN492   | NB basal medium                                                                     |         |         | 10L  |      |                                                |                |
| PMN613   | Nitsch and Nitsch basal salt mixture                                                |         |         | 10L  |      | Oxidising Agent/Irritant/Hygroscopic/Keep Cold |                |
| PMN616A  | Nitsch & Nitsch Basal Medium w/Vitamins 10L                                         |         | 10L     |      |      | Irritant                                       | Store at 2-6°C |
| PMN616B  | Nitsch & Nitsch Basal Medium w/Vitamins 50L                                         |         | 50L     |      |      | Irritant                                       | Store at 2-6°C |
| PMO139A  | Orchid maintenance/replate medium without charcoal and agar                         |         | 10 x 1L |      |      | Irritant                                       | Store at 2-8°C |
| PMO139B  | Orchid maintenance/replate medium without charcoal and agar                         |         | 10L     |      |      | Irritant                                       | Store at 2-8°C |
| PMO139C  | Orchid maintenance/replate medium without charcoal and agar                         |         | 50L     |      |      | Irritant                                       | Store at 2-8°C |
| PMO156A  | Orchid maintenance/replate medium with banana and charcoal, without agar            |         | 10L     |      |      | Irritant                                       | Store at 2-8°C |
| PMO156B  | Orchid maintenance/replate medium with banana and charcoal, without agar            |         | 10 x 1L |      |      | Irritant                                       | Store at 2-8°C |
| PMO156C  | Orchid maintenance/replate medium with banana and charcoal, without agar            |         | 50L     |      |      | Irritant                                       | Store at 2-8°C |
| PMO612A  | H1 Oat Medium                                                                       |         | 10 x 1L |      |      | Irritant                                       | Store at 2-8°C |
| PMO612B  | H1 Oat Medium                                                                       |         | 10L     |      |      | Irritant                                       | Store at 2-8°C |
| PMO753A  | Orchid multiplication medium with agar, without charcoal                            |         | 10 x 1L |      |      | Irritant                                       | Store at 2-8°C |
| PMO753B  | Orchid multiplication medium with agar, without charcoal                            |         | 10L     |      |      | Irritant                                       | Store at 2-8°C |
| PMO753C  | Orchid multiplication medium with agar, without charcoal                            |         | 50L     |      |      | Irritant                                       | Store at 2-8°C |
| PMP658A  | Orchid Maintenance Medium, Contains Charcoal and Agar                               |         | 10 x 1L |      |      | Irritant                                       | Store at 2-8°C |
| PMP658B  | Orchid Maintenance Medium, Contains Charcoal and Agar                               |         | 10L     |      |      | Irritant                                       | Store at 2-8°C |
| PMP658C  | Orchid Maintenance Medium, Contains Charcoal and Agar                               |         | 50L     |      |      | Irritant                                       | Store at 2-8°C |
| PMP668A  | Orchid maintenance medium with charcoal, without agar                               |         | 10 x 1L |      |      | Irritant                                       | Store at 2-8°C |
| PMP668B  | Orchid maintenance medium with charcoal, without agar                               |         | 10L     |      |      | Irritant                                       |                |

【お問い合わせ】 ナミキ商事株式会社 医薬化学品部

Tel: 050-5527-9830 (部署直通)

Tel: 03-3354-4065 (代表) Fax: 03-3354-3527

mail: info@namiki-s.co.jp <https://www.namiki-s.co.jp>



## Apollo バルク在庫品リスト

| カタログNo. | 品名                                                                         | Cas No. | Qty     | Qty2 | Qty3 | 備考                                                               | 保管条件           |
|---------|----------------------------------------------------------------------------|---------|---------|------|------|------------------------------------------------------------------|----------------|
| PMP668C | Orchid maintenance medium with charcoal, without agar                      |         | 50L     |      |      | Irritant                                                         |                |
| PMP723A | Orchid seed sowing medium                                                  |         | 10L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMP723B | Orchid seed sowing medium                                                  |         | 10 x 1L |      |      | Irritant                                                         | Store at 2-8°C |
| PMP723C | Orchid seed sowing medium                                                  |         | 50L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMP727  | Orchid seed sowing medium II                                               |         | 10L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMP748A | Orchid Maintenance / Replate Medium Contains Banana, Charcoal, and Agar    |         | 10L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMP748B | Orchid Maintenance / Replate Medium Contains Banana, Charcoal, and Agar    |         | 50L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMP748C | Orchid Maintenance / Replate Medium Contains Banana, Charcoal, and Agar    |         | 10 x 1L |      |      | Irritant                                                         | Store at 2-8°C |
| PMP782A | Phytotech orchid replate medium                                            |         | 10 x 1L |      |      | Irritant                                                         | Store at 2-8°C |
| PMP782B | Phytotech orchid replate medium                                            |         | 10L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMP782C | Phytotech orchid replate medium                                            |         | 50L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMP785A | Orchid replate medium                                                      |         | 10 x 1L |      |      | Irritant                                                         | Store at 2-8°C |
| PMP785B | Orchid replate medium                                                      |         | 10L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMP785C | Orchid replate medium                                                      |         | 50L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMP793A | Orchid multiplication medium without charcoal and agar                     |         | 10 x 1L |      |      | Irritant                                                         | Store at 2-8°C |
| PMP793B | Orchid multiplication medium without charcoal and agar                     |         | 10L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMP793C | Orchid multiplication medium without charcoal and agar                     |         | 50L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMQ673A | Quoirin and Lepoivre basal salt mixture 10 x 1L                            |         | 10 x 1L |      |      | Irritant                                                         | Store at 2-6°C |
| PMQ673B | Quoirin and Lepoivre basal salt mixture 10L                                |         | 10L     |      |      | Irritant                                                         | Store at 2-6°C |
| PMQ673C | Quoirin and Lepoivre basal salt mixture 50L                                |         | 50L     |      |      | Irritant                                                         | Store at 2-6°C |
| PMS816A | Schenk and Hildebrandt basal salt mixture                                  |         | 10 x 1L |      |      | Oxidising Agent/Irritant                                         | Store at 2-8°C |
| PMS816B | Schenk and Hildebrandt basal salt mixture                                  |         | 10L     |      |      | Oxidising Agent/Irritant/Hygroscopic/Store under Argon           | Store at 2-8°C |
| PMS816C | Schenk and Hildebrandt basal salt mixture                                  |         | 50L     |      |      | Oxidising Agent/Irritant/Hygroscopic/Store under Argon/Keep Cold | Store at 2-8°C |
| PMT839A | Terrestrial (Cypripedium) Orchid Medium                                    |         | 10 x 1L |      |      | Irritant                                                         | Store at 2-8°C |
| PMT839B | Terrestrial (Cypripedium) Orchid Medium                                    |         | 10L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMT839C | Terrestrial (Cypripedium) Orchid Medium                                    |         | 50L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMT842A | Terrestrial (Cypripedium) Orchid Medium, contains 600 mg/L Calcium nitrate |         | 10 x 1L |      |      | Irritant/Moisture Sensitive                                      | Store at 2-8°C |
| PMT842B | Terrestrial (Cypripedium) Orchid Medium, contains 600 mg/L Calcium nitrate |         | 10L     |      |      | Irritant/Moisture Sensitive                                      | Store at 2-8°C |
| PMT842C | Terrestrial (Cypripedium) Orchid Medium, contains 600 mg/L Calcium nitrate |         | 50L     |      |      | Irritant/Moisture Sensitive                                      | Store at 2-8°C |
| PMT849A | Terrestrial (Cypripedium) Orchid Medium, contains 400 mg/L Calcium nitrate |         | 10 x 1L |      |      | Irritant/Moisture Sensitive                                      | Store at 2-8°C |
| PMT849B | Terrestrial (Cypripedium) Orchid Medium, contains 400 mg/L Calcium nitrate |         | 10L     |      |      | Irritant/Moisture Sensitive                                      | Store at 2-8°C |
| PMT849C | Terrestrial (Cypripedium) Orchid Medium, contains 400 mg/L Calcium nitrate |         | 50L     |      |      | Irritant/Hygroscopic/Store under Argon                           | Store at 2-8°C |
| PMV505A | Vacin and Went modified orchid basal salt mixture                          |         | 10 x 1L |      |      | Irritant                                                         | Store at 2-8°C |
| PMV505B | Vacin and Went modified orchid basal salt mixture                          |         | 10L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMV505C | Vacin and Went modified orchid basal salt mixture                          |         | 50L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMV895A | Vacin and Went modified orchid medium with agar and sucrose                |         | 10 x 1L |      |      | Irritant                                                         | Store at 2-8°C |
| PMV895B | Vacin and Went modified orchid medium with agar and sucrose                |         | 10L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMV895C | Vacin and Went modified orchid medium with agar and sucrose                |         | 50L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMW898A | White's basal salt mixture                                                 |         | 10 x 1L |      |      | Irritant                                                         | Store at 2-8°C |
| PMW898B | White's basal salt mixture                                                 |         | 10L     |      |      | Irritant                                                         | Store at 2-8°C |
| PMW898C | White's basal salt mixture                                                 |         | 50L     |      |      | Irritant                                                         | Store at 2-8°C |

【お問い合わせ】 ナミキ商事株式会社 医薬化学品部

Tel: 050-5527-9830 (部署直通)

Tel: 03-3354-4065 (代表) Fax: 03-3354-3527

mail: info@namiki-s.co.jp <https://www.namiki-s.co.jp>

