


ACROBiosystem's PD1 protein collection featuring full-length PD1 protein

ACROBiosystems newly added the world's first full length human PD1 protein to an already extensive archive of PD1 family proteins. The full length protein demonstrates a much higher binding affinity against PDL1 than the PD1 ECD protein.



A representative ELISA result showed that immobilized full length human PD1 (PD1-HC214, 20 ng/well) can bind with human PDL1 with a linear range of 1.2-78ng, which is much higher than the PD1 ECD protein (linear range 250ng-1,250ng).

[Learn More](#)

We also provide a wide variety of PD1 family proteins from different species. This includes the exclusive PrecisionAvi series of AviTag™ biotinylated PD1, PDL1 and PDL2, the HPLC-verified products, and the high throughput PD1:PDL1 pathway inhibitor screening assay pair.

Molecule	Cat. #	Product Description
PD-1	PD1-HC214	Human PD-1 Protein, Tag Free, Full Length
PD-1	PD1-H82E4	Biotinylated Human PD-1 (recommended for biopanning), Avi Tag (Avitag™)
PD-1	PD1-H82F2	Biotinylated Human PD-1, His Tag & Fc Tag, Avi Tag (Avitag™)
PD-1	PD1-M8259	Biotinylated Mouse PD-1, Fc Tag, ultra sensitivity (primary amine labeling)
PD-1	PD1-M82F4	Biotinylated Mouse PD-1, Fc Tag, Avi Tag (Avitag™)
PD-L1	PD1-H82E5	Biotinylated Human PD-L1 (recommended for biopanning), Avi Tag (Avitag™)
PD-L1	PD1-H82F3	Biotinylated Human PD-L1, His Tag & Fc Tag, Avi Tag (Avitag™)
PD-L1	PD1-M82F5	Biotinylated Mouse PD-L1, Fc Tag, Avi Tag (Avitag™)
PD-L2	PD2-H82E8	Biotinylated Human PD-L2 (recommended for biopanning), Avi Tag (Avitag™)
PD-L2	PD2-H82F6	Biotinylated Human PD-L2, Fc Tag, Avi Tag (Avitag™)
PD-1	PD1-H5221	Human PD-1 (HPLC-verified)
PD-1	PD1-H5257	Human PD-1, Fc Tag (HPLC-verified)
PD-L1	PD1-H5229	Human PD-L1 (HPLC-verified)
PD-L1	PD1-H52A3	Human PD-L1, mouse IgG1 Fc Tag (HPLC-verified)

PD-L2	PD2-H5251	Human PD-L2, Fc Tag (HPLC-verified)
PD-1	PD1-H5284	Human PD-1, Strep Tag
PD-1	PD1-C5223	Cynomolgus PD-1
PD-1	PD1-C5254	Cynomolgus PD-1, Fc Tag
PD-1	PD1-M5228	Mouse PD-1
PD-1	PD1-M5259	Mouse PD-1, Fc Tag
PD-1	PD1-R5253	Rat PD-1, Fc Tag
PD-1	PD1-R52H2	Rat PD-1
PD-1	PD1-R52H0	Rabbit PD-1
PD-1	PD1-RB5251	Rabbit PD-1, Fc Tag
PD-1	PD1-H9221	C13 and N15-labeled, Human PD-1
PD-1	PD1-HF257	FITC-Labeled Human PD-1, Fc Tag
PD-L1	PD1-H5258	Human PD-L1, Fc Tag
PD-L1	PD1-H5282	Human PD-L1, Strep Tag
PD-L1	PD1-C5253	Cynomolgus PD-L1, Fc Tag
PD-L1	PD1-C52H4	Cynomolgus PD-L1
PD-L1	PD1-M5220	Mouse PD-L1
PD-L1	PD1-M5251	Mouse PD-L1, Fc Tag
PD-L1	PD1-M52A2	Mouse PD-L1, mouse IgG2a Fc Tag
PD-L1	PD1-H9229	C13 and N15-labeled, Human PD-L1
PD-L2	PD2-H5220	Human PD-L2
PD-L2	PD2-H52A5	Human PD-L2, mouse IgG1 Fc Tag
PD-L2	PD2-M5254	Mouse PD-L2, Fc Tag
PD-L2	PD2-M52E3	Mouse PD-L2
PD-1 & PD-L1	EP-101	PD-1 [Biotinylated] : PD-L1 Inhibitor Screening ELISA Assay Pair

STAY CONNECTED

