## ASINEX Signature Libraries

## **SL-10.** Purine-based Nucleosides

The enzyme co-factor SAM/SAH is found in all branches of life, widely exploited by protein methyltransferases for the transfer of a methyl group to various molecules of biological significance. A number of these enzymes are drug targets in the field of epigenetics [1].

There is also wide-ranging interest in SAM/SAHdependent enzymes beyond epigenetics, in the fields of antibacterials, anti-virals, anti-parasitics, oncology (which includes both epigenetic and non-epigenetic mechanisms), agrochemicals, and in industrial chemical processing.

Compounds targeting the SAM/SAH site of these enzymes offer great potential, yet very few such molecules are

currently described in literature or commercially available in compound collections.

To address this market need, we have used a modular approach to synthesize a library of purine-based nucleosides. The main challenge was to develop stereo- and enantio-selective methods for the synthesis of nucleotide-like cores which both mimic adenosine and are metabolically stable.



## Signature Library 10

Formats	Supplementary Information
80 compounds per plate	SL#10_PurineNucleoside_04-16.sdf
0.1 mg; 1 mg; 2 mg dry film/powder	
0.1 µmol; 1 µmol DMSO solutions	

asinex.com

## References:

1. ACS Chem Biol. 2012 Mar 16;7(3):443-63. doi: 10.1021/cb200519y.

<u>Contact us</u>: USA: Japan: Europe/Global:

+1 336 721 1617 +81-80-3401-9097 mparisi@asinex.com sota@asinex.com lsadovenko@asinex.com